

Checker®

Checker®

Checker®


Checker®

- Stir gently and wait until the display stabilizes.
- For best results, recalibrate periodically.
- NEVER IMMERS THE ELECTRODE UP TO THE CONNECTOR CLEAN AND DRY.
- After use, rinse the electrode with water to minimize contamination.
- Store the electrode with a few drops of **HI 70300 Storage Solution** in the protective cap.
- DO NOT USE DISTILLED OR DEIONIZED WATER FOR STORAGE PURPOSES.
- Always replace the protective cap after use.



OPERATION:

- Do not be alarmed if white crystals appear around the cap. This is normal with pH electrodes and they dissolve when rinsed with water.
- If the electrode is dry, soak it in tap water for a few minutes, prior to use.
- Connect the electrode to the meter.



- Switch the **Checker®** on.
- Remove the protective cap and immerse the tip of the electrode (bottom 4 cm/1½") into your sample.

Calibration Manual two points

Electrode: combination stick pH electrode

Checker® 1: HI 1270 (included)

Checker® 2: HI 1207 (included)

Checker® 3: HI 1208 (included)

Environment 0 to 50°C (32 to 122°F); 95% RH max.

Battery Type 2 x 1.5V alkaline approx. 3,000 hours of continuous use

Dimensions 66 x 50 x 25 mm (2.6 x 2 x 1")

Weight (meter) 70 g (2.5 oz.)


INITIAL PREPARATION:

The pH electrode is shipped dry. Prior to using the **Checker®**, remove the protective cap and condition the electrode by soaking the tip (bottom 4 cm/1½") in pH 7.01 buffer solution for a couple of hours. Then follow the calibration procedure below.

Checker® 1: with HI 1270 screw-type pH electrode

Checker® 2: with HI 1207 screw-type pH electrode

Checker® 3: with HI 1208 BNC-type electrode



SPECIFICATIONS:

Range 0.00 to 14.00 pH

Resolution 0.01 pH

Accuracy ±0.2 pH (@20°C/68°F)

Typical EMC Deviation ±0.1 pH